





AT LISA tri 839MP preloaded
Technical Specifications



We make it visible.

AT LISA tri 839MP preloaded

AT LISA® tri 839MP preloaded



Optic Design	Trifocal, diffractive, +3.33 D near add and +1.66 D intermediate add at the IOL plane, aspheric (aberration correcting)
Material	Hydrophilic acrylic (25 %) with hydrophobic surface
Optic Diameter	6.0 mm
Total Diameter	11.0 mm
Haptic Angulation	0°
Lens Design	Single-piece, MICS
Incision Size	1.8 mm
Company Labeled A-Constant ¹	118.6
Diopter Range	From 0.0 to +32.0 D 0.0 to +32.0 D, 0.5 D increments
ACD	5.32
Implantation in	Bag
Injector / Cartridge Set ²	BLUEMIXS 180
Indications	Presbyopia correction in patients with or without cataract (Prelex or clear lens exchange)

¹Please refer to our web pages for optimized A-Constants. ²Please refer to our web pages for the most up-to-date references.



0297 AT LISA tri 839MP



0459 BLUEMIXS 180 injector



Carl Zeiss Meditec AG
Goeschwitzer Str. 51–52
07745 Jena
Germany
www.meditec.zeiss.com/lisa-tri
www.meditec.zeiss.com/contacts



Carl Zeiss Meditec SAS
Avenue Paul Langevin, BP5
17053 La Rochelle Cedex 9
France
www.meditec.zeiss.com/iol
www.meditec.zeiss.com/contacts

ID/01167 02/12 GB

AT LISA and BLUEMIXS are registered trademarks of Carl Zeiss Meditec. The contents of the brochure may differ from the current status of approval of the product in your country. Please contact our regional representative for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. Printed on elemental chlorine-free bleached paper.
© 2012 by Carl Zeiss Meditec AG. All copyrights reserved.



000000-1982-735



02/2012

